



ACCESSORIES FOR HOSPITAL BED

ATTENDANT BED

Code : KMMEI- 123

Description

- Frame made of rectangular CRCA pipes.
- Bed top made of CRCA perforated sheet.
- Frame is mounted on PVC stumps.
- The complete metal components are pre treated and epoxy powder coated.
- Frame Size : 1870L x 600W x 450H
- Note : Mattress optional



BABY CRADLE ON STAND

Code : KMMEI- 124

Description

- Frame made of M,S. Tubular
- Legs with rubber stumps
- Detachable Baby cradle
- Epoxy Powder Coated



FOLDABLE BACKREST

Code : KMMEI- 125 B

Description

- Power Coated Backrest
- Size : 605x685x605/410mm
- Packing size : 700x630x100mm (1pcs)
- N.W.:3KG



BED SHIFTING PLATE

(Move Care Plat)

Code : KMMEI- 125 S

Description

- Structure: Foldable
- Material: Polypropylene
- Size(mm): Expansion Length(620-660mm), Width(340-360mm), Thickness \leq 4.5mm.
- 4.Move Care Plate surface has a certain sliding friction coefficient, The underside of material has stop slippery characteristics and equipped with the slide feature.
- To ensure the products strength and safety performance. Between Move care plate and transfer objects, when the distance of the contact point have reach 200 mm, the bearing strength is still not less than 100kgs.
- Easy to carry, the folded width is not more than 130 mm.
- The contact surface have no exposed cutting edge, cutting edges, the blade and hairy spine. 8.The injection mold with smooth surface, without defect, without cave, uniform color.



MOSQUITO NET POLES

Code : KMMEI- 126

Description

- Made of M.S. Tubular Steel.
- Epoxy Powder Coated



ACCESSORIES FOR HOSPITAL BED

HEAD & FOOT BOARD (ABS)

Code : KMMEI- 127A



ADJUSTING FOOD TABLE FOR BED (ABS)

Code : KMMEI- 127B

COLLAPSIBLE SAFETY SIDE RAILING (PAIR)

Code : KMMEI- 128

Description

- Size : 1500 / 1300 mm
- Aluminium Made



DUCK TYPE SAFETY SIDE RAILING

Code : KMMEI- 128A



DROP SIDE RAILING (PAIR)

Code : KMMEI- 128B

Description

- Made of MS Tubular
- Epoxy powder Coated





ACCESSORIES FOR HOSPITAL BED

I.V ROD S.S.

Code : KMMEI- 128 C



PATIENT CASE SHEET HOLDER

Code : KMMEI- 129 (Double Flap)

Description

- Double Flap
- Made of Aluminium Sheet



2" WHEEL CASTORS

Code : KMMEI- 130A

- Luxurious Silent Double-Sided Castor



6" WHEEL CASTORS

Code : KMMEI- 130C

- Luxurious Silent Double-Sided Castor



4" WHEEL CASTORS

Code : KMMEI- 130B

- Luxurious Silent Double-Sided Castor