



## CHEST FREEZER & ICELINE REFRIGERATOR

## ICE-LINED REFRIGERATOR

**Code : KMMEI- 650**

### Description

#### Patented Technology

- Heated pipe for better temperature uniformity

#### CFC-free

- CFC-free high-density foam insulation

#### Temperature Control

- Equipped with digital solar powered temperature display to measure controlled inside temperature range from 2°C– 8°C

#### Ergonomic Design

- Safety lock to avoid unauthorized access control, equipped with handles on the sides



HBCD-90



Solar Energy Display Panel



HBC-80



Solar Energy Display Panel



### Specifications

Model	HBC-80	HBC-150	HBC-260	HBCD-90	HBC-120	HBC-240
WHO PQS Code	Refrigerator	Refrigerator	Refrigerator	Refrigerator + Freezer	Refrigerator	Refrigerator
Cabinet Type	E003/089	E003/088	E003/087	E003/097	Coming soon	Coming soon
Ambient	Chest	Chest	Chest	Chest	Chest	Chest
Temperature (oC)	5~43	5~43	5~43	5~43	5~43	5~43
Gross Volume (L)	80	150	260	Refrig.: 58 Freezer: 32	120	240
Vaccine Storage Capacity	61	122	211	37,5	100	200
Exterior Dimension	788*717*872	1128*717*872	1647*717*940	1128*720*875	865*825*1422	865*825*1815
Holdover Time at 43oC	59hrs 58mins	60hrs 50mins	More than 60hrs	137hrs 47mins	100hrs	87hrs
Noise Level (dB(A))	<40	<40	<40	<40	<40	<40
Freeze Protection Level	A	A	A	A	A	A
Optional	Temperature Logger Automatic Voltage Stabilizer Remote Temperature Monitoring Device					