

DIGITAL THERMOMETER, MUAC TAPE

DIGITAL THERMOMETER

Description

Code : KMMEI- 327

- Rigid tip
- Last reading recall
- Dual scale, (DF) Fahrenheit and (°C) Celsius both, Beeps,
- Fever alarm in 60 Seconds
- Auto-off



DIGITAL THERMOMETER

Description

Code : KMMEI- 328

- Fever Line indicator
- Jumbo LCD,
- Fast read in 20 seconds
- Flexible Tip, Waterproof
- Last reading recall
- Dual scale, (0F) Fahrenheit and (0C) Celsius both,
- Break resistance,
- Beeps
- Fever alarm
- Auto-off



With outer hard cover case



DIGITAL BLOOD GLUCOSE MONITOR

Description

Code : KMMEI- 335

- BG-202 Blood Glucose Monitoring System
- Assay method Glucose oxidase biosensor
- Measurement range 2.2-33.3 mmol/L (40-600 mg/dL)
- Sample Fresh capillary whole blood
- Unit of measure; mmol/L or mg/dL; switchable
- Memory ; 360 blood glucose results with date and time
- LCD screen size: 34.6\*31.1 mm (L\*H)
- Sample size ; Approximate 1 microlitre
- Response time ; 5 seconds
- Coding ; No coding
- Battery ; One 3.0V CR2032 lithium battery, replaceable
- Battery life ; Approximately 1,000 tests
- Standard: ISO 15197. & EN 61326,
- Calibration Plasma equivalent glucose values
- Automatic power off: Approx. 3 minutes after last user action



Description

Code : KMMEI- 337

MUAC TAPE FOR PREGNANT WOMEN

- The AnthroFlex mid-upper arm circumference (MUAC) tape measure is an economical solutions that allows you to quickly and easily assess nutritional risk of adults over 18 years, and postpartum and pregnant women.
- ASSESS MALNUTRITION RISK quick and easily in adults over 18 years, and postpartum and pregnant women, based on mid-upper arm circumference (MUAC) reading. EASY TO USE Red represents severe malnourishment.
- Yellow represents at-risk for malnourishment, and Green represents proper nourishment.



MUAC TAPE FOR CHILD ARM

Description

Code : KMMEI- 338

MUAC measuring tape is suitable for measuring child's Middle Upper Arm Circumference, (MUAC) with range up to 26.5 cm. Graduated with 1 mm precision with thicker line at 21.0 cm. Accuracy: ± 1 mm of the maximum measurement (26.5 cm), Front side: Color-coded as follows: Red : from 0 to 11.5 cm, Yellow from 11.5 to 12.5 cm, Green from 12.5 to 26.5 cm. Material: Non-tear stretch-resistant polyethylene material. Minimum thickness 0.3mm,



MUAC TAPE FOR ADULTS

Description

Code : KMMEI- 339

Adults' Mid Upper Arm Circumference (MUAC) measuring tape with a range of up to 56cm. Suitable for measuring adult's Middle Upper Arm Circumference (MUAC, range to 56 cm). Graduated with 1 mm precision. Accuracy: ± 1 mm of the maximum measurement (50 cm), Material: Non-tear stretch-resistant polyethylene material. Minimum thickness 0.3mm.

