



FACIAL STEAMER, DIGITAL BP

FACIAL STEAMER

Code : KMMEI- 322-SP

Description

- 220-230V,130W
- Accessories: face mask, nose mask, measuring cup, manual,
- Using time: around 10 ~ 15 minute,
- Material: face/nose mask-PVC, main body PP, pot-stainless iron,
- 2 meter non heating, good quality wire,
- EU plug



ARM-TYPE DIGITAL BLOOD PRESSURE MONITOR

Code : KMMEI- 325



Leading MDI Technology
Faster measurement
Upgraded experience



Description

- LCD size 66*43mm
- Blood pressure result indicator
- Low Battery Detection
- Automatic off
- Power Source: 3xAAA battery or Type -C cable(not included)
- Standard 22-36cm cuff
- Drawstring bag
- Back light
- Talking
- Bluetooth



INFRARED THERMOMETER

Code : KMMEI- 326

Description

- Measure on forehead, Safe and Accurate,
- No probe covers needed
- 10 reading memory
- 3 second reading
- Dual scale, (°F) Fahrenheit and (°C) Celsius both,
- Compact ergonomic size
- Beeps when ready
- Jumbo display with backlight
- Replaceable battery
- Auto-off
- Response Time: 3 seconds
- Fever Alarm: over 37.80C (100.40F)
- Unit Weight: Approx. 84.0 grams
- Range: 34.00C- 43.00C(89.6 0F-109.4 0F)
- Unit Dimension: 15.4X4.1X4.5cm
- Battery: DC3V(2xAAA battery)

