

GENERAL EQUIPMENT

SYRINGE CUTTER & NEEDLE DESTROYER

Code : KMMEI- 257

Description

- It can destroy needle and Syringe (2ml to 10ml Capacity)
- Power Supply 220V + 10% 50Hz AC



SYRINGE & NEEDLE CUTTER MANUAL

Code : KMMEI- 257 M

Description

No-Sharps (NS-800DX) Disposable Container of 800 ml, Cuts The syringe hub & Needle in single action, Has Stainless steel frame & facility for wall mounting (100,000 cuts Life)

PLASTIC DUSTBIN

Code : KMMEI- 258

Description

- Available in Many sizes
- Material : PP/ CP Injection Moulded
- Dimension : L 26.8 x W 24.2 x H 28.5 CM
L 30.5 x W 30.5 x H 35.6 CM
L 45.5 x W 29.5 x H 49.5 CM
- Available in colour coding options of Red, Blue, Black, Yellow & Green



DUST BIN / WASTE BIN

Code : KMMEI- 259

Description

- Available in Stainless Steel, and PVC,
- Available in different sizes and shapes,



GENERAL EQUIPMENT

MICRO FOGGER

Code : KMMEI- 260

Description

- Max. area coverage as compared to any other conventional fogger.
- 304 Grade Stainless Steel Body.
- Chemical Level Indicator.
- Comes with pre-fitted liquid flow control valve for ease of use.
- 100% Convenient, Portable & Easy to use.
- Body & Tank : 304 Grade Stainless Steel
- Tank Capacity : 5 Litre of Chemical Storage
- Blower Motor : >26000 RPM Motor
- Nozzle : 1 Counter Rotation Vortex Nozzle
- Chemicals : Can atomize both oil-based and water-based liquids
- Control Valve : One turn. Brass Valve
- Particle Size : 5-30 μ , adjustable. Liquid viscosity and density can affect particle size
- Range: Visible fog, 20-30 ft. (2.5m)
- Input: 230V, 50Hz, Single Phase, AC Supply
- Dimensions: [33 x 41 x 23.5 cm] L x H x Dia
- Net. Weight: 6.50 Kg (Approx.)



AEROSOL O.T. DISINFECTOR

Code : KMMEI- 261

Description

- Made of SS
- Capacity 5Ltr.
- OT Fumigator
- For the purpose of disinfecting or
- eradicating pests in Hospital OT.



OXYGEN HOOD JOINTLESS

Code : KMMEI- 262

Description

- Single piece round shape up of autoclave polycarbonate material
- Trauma free silicone neck, adjustable flap.
- Bilateral oxygen nozzle, prevents direct flow of cold oxygen on patient head.
- Height app. : 7.5 inches
- Width app.: 8.5 inches

GENERAL EQUIPMENT

VACUUM EXTRACTOR GLASS TYPE

Code : KMMEI- 265

Description

- 2 silicon cups, 50mm, & 60mm with metal traction handle,
- Glass Bottle with Stand,
- 70mm Silicon cup is also available on request,



VACUUM EXTRACTOR POLYCARBONATE TYPE

Code : KMMEI- 266

Description

- 2 silicon cups, 50mm, & 60mm with metal traction handle,
- Autoclavable polycarbonate jar with pressure gauze,
- 70mm Silicon cup is also available on request,



SILICON VACUUM CUP WITH METAL HANDLE

Code : KMMEI- 267



The Ayre spatula is a device used to collect Pap smear. It is a wooden spatula with U shaped openings on one side and a flat surface on another. The broad end is for vaginal sample collection and the narrow end is for cervical sample collection.



AYRES SPATULA

Code : KMMEI- 270

Description

- Ayres Spatula with Cytobrush,
- Sterilized packing,

AIR RING/ SEAT CUSHION

Code : KMMEI- 272

Description

- An inflatable rubber ring used as a cushion' to relieve pressure on prominent bony points (as in helping to prevent bedsores),
- Available in Different Sizes.



AUTOCLAVE STEAM STERILIZER TAPE

Code : KMMEI- 275

Description

- Autoclave Tape for steam sterilization,
- Size 19mm x 50 Mtr.
- Changes color after exposure to high temperature, typically 1210C in a steam autoclave,
- Technical criteria as per ISO standard
- Individual Packing,

