

HOSPITAL HOLLOWWARE & UTENSILS

**INSTRUMENT TRAY WITH COVER S.S.**

Code : KMMEI- 290

**Description**

- 9" x 6"
- 11" x 7"
- 10" x 8"
- 12" x 10"
- 15" x 12"
- 18" x 12"



**KIDNEY TRAY S.S.**

Code : KMMEI- 292

**Description**

- 6", 8", 10" & 12"
- Also available with cover



**DRESSING DRUM S.S. HEAVY JOINTED**

Code : KMMEI- 294

**Description**

- 9" x 9"
- 11" x 9"
- 11" x 5"
- 15" x 12"
- 11" x 18"



**DRESSING DRUM S.S. SEAMLESS**

Code : KMMEI- 296

**Description**

- 9" x 5"
- 9" x 9"
- 11" x 9"
- 11" x 5"
- 12" x 10"
- 14" x 5"
- 15" x 12"



**UV STERILIZER**

Code : KMMEI- 298

**Description**

- Tabletop Ultraviolet sterilization/disinfection device
- Available in different sizes,



**SURGICAL WASH BASINS**

Code : KMMEI- 300

**Description**

- 23cm, 28 cm, 32cm, 40cm
- Available in different sizes,





**HOSPITAL HOLLOWWARE & UTENSILS**

**BED PAN MALE/FEMALE**

**Code : KMMEI- 302**

**Description**

- Overall size: 375L x 310W x 70H mm



**BED PAN MALE/FEMALE**

**Code : KMMEI- 303**

**Description**

- Overall size: 350L x 285W x 97H mm



**BED PAN MALE**

**Code : KMMEI- 305**

**Description**

- Overall size: 405L x 285W x 75H mm



**BED PAN MALE / FEMALE WITH COVER**

**Code : KMMEI- 306**

**Description**

- Overall size: 350L x 285W x 97H mm



**URINAL MALE**

**Code : KMMEI- 308**

**Description**

- Capacity : 1030 ml



**URINAL FEMALE**

**Code : KMMEI- 312**

**Description**

- Capacity : 1010 ml



**SPITTING MUG**

**Code : KMMEI- 313**

**Description**

- Capacity : 450 ml
- Size : 100 x 50 mm



**FETOSCOPE**

**Code : KMMEI- 314**

**Description**

- Pinard Horn Fetoscope
- Used for listening to the heartbeat of a fetus.
- Available in Aluminium, Plastic, and Wooden,



**ALUMINIUM FETOSCOPE**

**Code : KMMEI- 314A**



**PLASTIC FETOSCOPE**

**Code : KMMEI- 314P**



**WOODEN FETOSCOPE**

**Code : KMMEI- 314W**



**HOSPITAL HOLLOWWARE & UTENSILS**

**PLAIN STERILIZER**

**Description**

- Are electrically fitted with 1Kw. -2Kw. single phase element of 220V/50HZ. as per sizes & three core wire.
- Are available with thermostat to control temperature of water from size 12ironwards.
- Are available with thermostat and tray lifting arrangement from size 12" onwards.
- Are available with auto-cut off element for all sizes.
- All sizes are approximate.
- Product can be available in various sizes as per customer requirement.

**Code : KMMEI- 315**



**Models/ Variants**

Sizes in Inches	Sizes in mm
8 X 4 X 3.5	205X100X90
8 X 5 X 4.5	205X125X120
10 X 5 X 4.5	255X125X120
12 X 6 X 6	305X155X150
14 X 6 X 6	355X155X150
16 X 6 X 6	405X155X150
16 X 8 X 8	405X205X200
18 X 8 X 8	460X205 X200
20 X 8 X 8	510X205X200
24 X 8 X 8	610X205X200
24 X 10 X 8	610X255X200
24 X 12 X 8	610X305X200
24 X 12 X 12	610X305X305



**Models/ Variants**

Sizes in Inches	Sizes in mm
12 X 6 X 5	300X155X130
16 X 6 X 5	410X155X130
17 X 8 X 6	430X205X150
20 X B X 6	510X205X150

**INSTRUMENT STERILIZER**

**Code : KMMEI- 317**

**Description**

- Product is completely fitted with thermostats & ray lifting arrangement but cane made available without thermostat & tray lifting arrangement.
- Product also available with plain immersion element or auto cut off. Product electrically fitted with 1Kw- 2kw. single phase element of 220V/50Hz. as per sizes & three core wire
- lifting arrangements can be made as per jointed instrument sterilizer.

**STERELIZER ELECTRIC S.S. AUTO CUT OFF**

**Description**

**Sizes :**

- 10" x 5" x 3"
- 12" x 6" x 4"
- 16" X 6" X 6"
- 18" X 8" X 6"

**Code : KMMEI- 318**



**MAGNIFYING GLASS WITH LED LIGHT**

**Code : KMMEI- 320**



**STEAMER SPECIAL TYPE**

**Code : KMMEI- 322-L**



**STEAMER MYST TYPE**

**Code : KMMEI- 322-M**

