



MEDICAL EQUIPMENT & ACCESSORIES

RING PASSARY SILICON DIFFERENT

Code : KMMEI- 384

Description

- Silicone Ring Pessary is very common for a first to second degree prolapsed uterus for urine incontinence. It also is used on an accompanying cystocele.
- Features:
Insertion is eased with the folding action of the ring, Soft and pliable, Easy to insert and remove, Material 100% Medical Grade Silicone, Available in Different Size,



MVA SET WITH 8 CANNULA & BOX

Code : KMMEI- 386

Description

We are offering a gamut of MVA Syringe having rate of infection as low as 0.5%. These products are used to induced abortion, which is a therapeutic procedure used after miscarriage. Offered products are sourced from renowned manufacturers ensuring high quality products. The range of products is highly efficient and works flawlessly. Patrons can avail these products from us in different sizes ranging 4 mm to 12 mm as per their needs and requirements.

Features: User friendly, Cost-effective, Reliable,

- Used as a method of induced abortion, a therapeutic procedure used after miscarriage.
- The rate of infection is lower than any other surgical abortion procedure at 0.5%.
- Cannulas are available in different sizes ranging from 4 mm-12 mm.

NEBULIZER BIG TYPE

Code : KMMEI- 389

Description

Advantages:

- Compressor with 3C, CE and TUV certifications;
- Noiseless device with high nebulizing efficiency;
- Small nebulizing particles (diameter: 1-5 micron for more than 60% of the particles);
- Portable and convenient design.

Type of protection

Against electric shock :Class II Equipment / Internally powered equipment

Degree of protection

Against electric shock: Type BF equipment
 Electric Power: ~110/220 V, 50/60 Hz;
 Power Consumption: ≤160 VA;
 Noise: ≤60 dB(A);
 Average Nebulizing Rate: ≥0.30 mL/min;
 Output Pressure: 60 ~130 kPa;
 Working Temperature/Humidity: 5~40°C/ 10~90% RH
 Size(mm) 285*180*110
 N.W(kg) 1.6
 Standard Accessories: Nebulizer (1)
 Air Tube (1)
 Angled Mouth Piece (1),
 Air Filter (3),
 Standard Adult Mask(1)



NEBULIZER PORTABLE TYPE

Code : KMMEI- 390

Description

- Compressor Nebulizer
- Noiseless device with high nebulizing efficiency, designed against electric shock, Type BF equipment
- small nebulizing particles (diameter: 1-5 micron for more than 60% of the particles,
- Portable and convenient (Teddy Bear/ Compact) design,
- Power consumption: 130VA+15%,
- Noise: <65dB (A)
- Average Nebulizing rate: >0.12 mL/min,
- Output Pressure: 60~130 kPa,



MEDICAL EQUIPMENT & ACCESSORIES

HOT-COLD PACK

Code : KMMEI-392

Description

- Cool pack help in relief pain, sports injuries, swelling and sore muscles, itches and sprains. It can use for headache or neck aches and pain. It helps to slow the nerve impulses in the area, which interrupts the pain-spasm reaction between the nerves.
- Heat can increase blood flow and help restore movement to injured tissue. Warmth can also reduce joint stiffness, pain, and muscle spasms. As with cold packs, heat packs have a role in easing pain from both acute and chronic injuries, such as sprains, strains, muscle spasms, whiplash, and arthritis.



OXYGEN FLOW METER WITH SPANNER

Code : KMMEI- 393

Description

- Gas Type : Medical Oxygen, Medical Air use
- Flow range 0-15 LPM
- Back pressure compensated flowmeters
- 3000 psi maximum inlet pressure
- Easy to read flow rate window
- Complies to CGA standards

INSULIN TRAVEL BAG

Code : KMMEI-395

Description

- BIG STORAGE - Stores two medication vials in the ice brick and the mesh pocket is able to hold other diabetic supplies as needed.
- COOLING COMFORT - The Micro Cooler allows up to 12 hours of cooling time so that you can relax and not worry about your insulin, medication, or other prescriptions that need to be cooled.
- DIABETIC NEEDS - Mesh elastic pocket for holding smaller supplies
- ICE BRICK - The ice brick is re-freezable so that you can find comfort again and again giving you more time to think about other priorities,
- TRAVEL SIZE - The product dimension are: 8 inches by 2 inches by 4 inches which allows the cooler to be stored in suitcases, backpacks, messenger bags, and more. It can be used for long car rides, over sea travel, and many other vacation spots.





MEDICAL EQUIPMENT & ACCESSORIES

VACCINE CARRIER DIFFERENT SIZE

Code : KMMEI- 396

Description

WHO-PQS Certified, Available in 3 sizes:

- Vaccine Storage Capacity: 0.9 Litres, Cold Life without opening, tested by: 21.5 hours at 43°C
- Vaccine Storage Capacity: 2.90 Litres, Cold Life without opening, tested by: 51.08 hours at 43°C
- Vaccine Storage Capacity: 1.67 Litres, Cold Life without opening, tested by: 35.7 hours at 43°C



COLD BOX DIFFERENT SIZE

Code : KMMEI- 397

Description

- WHO-PQS Certified, Available in 5 sizes:
- Vaccine Storage Capacity: 22.5 Litres, Cold Life without opening, tested by: 128.16 hours at 43°C
- Vaccine Storage Capacity: 22.45 Litres, Cold Life without opening, tested by: 140 hours at 43°C
- Vaccine Storage Capacity: 22.40 Litres, Cold Life without opening, tested by: 153.5 hours at 43°C
- Vaccine Storage Capacity: 8.0 Litres, Cold Life without opening, tested by: 84.07 hours at 43°C
- Vaccine Storage Capacity: 5.50 Litres, Cold Life without opening, tested by: 107.50 hours at 43°C



PEN TORCH

Code : KMMEI- 398

Description

- Powered by 2 AAA batteries (not included)
- Suitable for Indoor and outdoor use.
- Dimension:13.5cm(L) *1.2cm(Dia)(approx)
- Made of anodized aluminum
- Clip to hold it on your uniform, easy to carry



TONGUE DEPRESSOR

Code : KMMEI- 401

Description

- Wooden Tongue Depressors Non Sterile Non Splinter Chemical Free Wooden 6 inch High Grade Natural Birch.
- It is designed for strength and durability, made of 100% high quality wood and uniform in color, sturdy and constructed to be the smoothest tongue depressor in the mark.
- Used in medical practice to depress the tongue to allow for examination of the mouth and throat; for spreading ointments, mixing medications and even finger splints or in any beauty saloon as spatulas.
- These special Birch wood products are clean and smooth and has no odor or taste.



MEDICAL EQUIPMENT & ACCESSORIES

RUBBER HAND GLOVES

Code : KMMEI- 402

Description

The thickness of utility work gloves and its long cuffs provide excellent protection for all general cleaning and industrial tasks. Users get outstanding comfort from its silk lining. These rubber gloves have heat, chemical and oil resistance properties. Household Gloves, Medical & Surgical, Disposable Gloves, industrial, utility



RUBBER GUM BOOT

Code : KMMEI- 403

Description

- A High Quality Rubber made Gum Boot. Best use for Slippery and Muddy Working Area.
- High Quality Rubber
- Safety for Mud and Dirt
- Safety for Slippery Road
- Covers Toe to below Knee
- Thick and Sticking Sole for Slippery Working Area

McINTOSH RUBBER SHEET

Code : KMMEI- 404

Description

- Made of High Quality Natural Rubber
- Double faced colour Sheet Green/Red or Red/Blue
- Soft & High quality Breaking & Tearing Strength
- Used in Hospitals, Homes for lying under Babies on the bed & General use.
- One Roll of 10 meters.
- Color : Red/Blue or Green/Red (Dispatch randomly as per available color)



SURGICAL BLADE

Code : KMMEI- 405

Description

- Sterilized Surgical Blade Different Size:
- # 10, 11, 12, 15, 18, 20, 21, 22, 23, 24,
- Made of Carbon Steel, also available in Stainless Steel.
- Box of 100 pieces.

SURGICAL SCALPE HANDLE

Code : KMMEI- 409

Description

- SUPERIOR STAINLESS-STEEL/CARBON STEEL BLADE: Equipped with tempered Stainless-Steel blades with a plastic handle, disposable scalpel knives allow for heavy duty use and a comfy grip for easy control.
- INDIVIDUALLY WRAPPED & STERILE: Works amazingly for dermaplaning or small aesthetic procedures such as callus removing, dead skin-skin tag removal and scraping, acne removal and more.
- HIGHLY VERSATILE: Useful for premed student supplies, dissections at school/ university, taxidermy, clay modeling or any arts crafts project that requires a sharp scalpel.
- PLASTIC HANDLE WITH MEASUREMENT: The handle is designed to offer a good grip and allows for delicate control when required.



MEDICAL EQUIPMENT & ACCESSORIES

BP MACHINE ANEROID ADULT

Code : KMMEI- 425

Description

- Aneroid Sphygmomanometer Blood Pressure Cuff for Adult
- Features a standard inflation bulb, standard air release valve and calibrated cuff closure, to insure a secure fit and accurate reading.
- Applies the mechanical manometer to measure the blood pressure
- Easy to maintain, Great for personal home use, travel and office.
- Works equally well in hospitals, nursing homes or home health care,
- 20~300mmHg manometer attains the accuracy of +/- 3 mmHg with zero pin stop,
- Easy to Read an identifiable high-contrast dial, easy to read the numbers.
- The gauge dial is outfitted with a heavy-duty convenient clip, so it can be easily attached to the gauge holder on the cuff or clothing.
- The inflation bulb is crack-resistant, non-sticking, high-density, sturdy and durable, The easy-to-use large screw-type valve delivers precise deflation control rate.
- Portable Blood Pressure Case/Bag is included for safe and easy storage at home and on the go,




BP MACHINE ANEROID CHILD

Code : KMMEI- 426

Description

- Aneroid Sphygmomanometer Blood Pressure Cuff for Child.
- Features a standard inflation bulb, standard air release valve and calibrated cuff closure, to insure a secure fit and accurate reading.
- Applies the mechanical manometer to measure the blood pressure
- Easy to maintain, Great for personal home use, travel and office.
- Works equally well in hospitals, nursing homes or home health care,
- 20~300mmHg manometer attains the accuracy of +/- 3 mmHg with zero pin stop,
- Easy to Read an identifiable high-contrast dial, easy to read the numbers.
- The gauge dial is outfitted with a heavy-duty convenient clip, so it can be easily attached to the gauge holder on the cuff or clothing.
- The inflation bulb is crack-resistant, non-sticking, high-density, sturdy and durable, The easy-to-use large screw-type valve delivers precise deflation control rate.
- Portable Blood Pressure Case/Bag is included for safe and easy storage at home and on the go,



BP MACHINE ANEROID INFANT

Code : KMMEI- 427

Description

- Aneroid Sphygmomanometer Blood Pressure Cuff for Infant.
- Features a standard inflation bulb, standard air release valve and calibrated cuff closure, to insure a secure fit and accurate reading.
- Applies the mechanical manometer to measure the blood pressure
- Easy to maintain, Great for personal home use, travel and office.
- Works equally well in hospitals, nursing homes or home health care,
- 20~300mmHg manometer attains the accuracy of +/- 3 mmHg with zero pin stop,
- Easy to Read an identifiable high-contrast dial, easy to read the numbers.
- The gauge dial is outfitted with a heavy-duty convenient clip, so it can be easily attached to the gauge holder on the cuff or clothing.
- The inflation bulb is crack-resistant, non-sticking, high-density, sturdy and durable, The easy-to-use large screw-type valve delivers precise deflation control rate.
- Portable Blood Pressure Case/Bag is included for safe and easy storage at home and on the go,



MEDICAL EQUIPMENT & ACCESSORIES

STETHOSCOPE ADULT

Code : KMMEI- 428

Description

- Premium Stethoscope Excellent acoustic performance, the dual head stethoscope is ideal for nurses, medical students and professional doctors, helping to identify heart, lung and other body sounds.
- IDEAL FOR DIAGNOSTICS - An excellent general-purpose stethoscope for all manner of medical professionals. Designed to assist doctors, nurses, and EMTs with patient diagnostics and evaluations.
- DUAL-HEAD STETHOSCOPE - Designed for patients of all ages, the stethoscope is made for nurses, doctors and medical students looking for a high-quality stethoscope. The eartubes are both flexible and adjustable to ensure the most comfortable fit for the user. The solidly constructed chest piece allows for superior contact with the patient to provide more accurate readings and evaluations.
- ANATOMICALLY DESIGNED - The specially designed headset is made to match the angle of the ear canal for a superior fit that allows for better performance. Large 27-inch length allows you to easily accommodate all patients from infants to adults.
- COSPLAY AND ROLE PLAY - Heighten the fun and realism of your role play and cosplay activities by adding this functioning stethoscope to your outfit to take actual readings. It is durable and available in multiple colors to easily match your clothing.



STETHOSCOPE CHILD

Code : KMMEI- 429



STETHOSCOPE INFANT

Code : KMMEI- 430



STETHOSCOPE STAINLESS STEEL

Code : KMMEI- 431



KANGAROO MOTHER CARE (KMC) CHAIR

Code : KMMEI - 432

Kangaroo mother care (KMC) is a method of care of preterm or low birth weight (LBW) infants by placing them in skin to skin (STS) contact with mother or other caregiver in order to ensure optimum growth and development of the infant. Initially devised as an alternative to conventional technology-based care, KMC is now considered as a standard of care for LBW infants for all settings as an adjunct or alternative to conventional technology based care.

Metal Frame, AM Epoxy Powder Coated, Foldable, Waterproof Mattress, Mattress Color Blue, Black, or other options are available, Height Adjustable, Reclining position, adjustable hand rest, with leg support,





MEDICAL EQUIPMENT & ACCESSORIES

ARI SOUND TIMER

Code : KMMEI- 435

Description

- Respiratory Rate Timer
- The Respiratory Rate Timer has been carefully made of provide a life span of 10,000 counting cycles of 60 seconds, under normal usage.
- Specifications:
- Timing duration : 60 Seconds
- Timing accuracy : Within 0.1 seconds
- Timebase : Crystal Oscillator OF 32758 hZ
- Neck Cord : Approximately one meter long
- Operating temperature : 0 to 50 degrees C
- Storage temperature : -10 to 50 degrees C



X-RAY VIEW BOX

Code : KMMEI- 436

Description

This innovative x-ray illuminator finds a vast application in hospitals and nursing homes. Ideal for viewing x-ray films and MRI, this offered illuminator is widely acclaimed in the market. Known for its long life span, it is immensely appreciated. Further, to suit the convenience of the customers,

Features:

- No UV emission, Initial brightness can be set
- Reliable magnetic nipper,

Specifications:

- SINGLE SCREEN, DOUBLE SCREEN, & FOUR SCREEN,
- Power Consumption: 60W ~ 100W
- Light Source: LED
- Power Supply Input
- AC100-240V, 50/60Hz
- Brightness (LUX) >5,500,



URIN BAG

Code : KMMEI- 438

Description

- Made of medical grade material, it is soft, flexible, non-toxic and safe.
- Clear bag with detailed scale printed on. Every user could know the exact volume of waste liquid
- With aid part, even patient can finish the task easily and quickly.
- It is not embarrassing to collect or discharge urine for patients.



PEAK FLOW METER

Code : KMMEI- 440

Description

Peak Flow Meter is a portable, easy-to-use device that measures how well your lungs are able to expel air. By blowing a quick blast of air through a mouthpiece on one end, the peak flow meter can measure the force of air in liters per minute and give you a reading on a built-in numbered scale.



MEDICAL EQUIPMENT & ACCESSORIES

**PEAK FLOW METER
DISPOSABLE MOUTHPIECES**

Code : KMMEI- 441

Description

Cardboard disposable mouthpieces for use with peak flow meters, large spirometers and breathing training devices. These wax coated, universal disposable mouthpieces have a 30mm diameter. As they are disposable, this helps to maintain hygiene standards, and eliminates cross-infection issues.



ROOM THERMOMETER DIGITAL

Description

Code : KMMEI- 445

This indoor thermometer is a must have to maintain right humidity and temperature level for sensitive humans, animals, plants and more. As you know, too dry or too wet, too hot or too cold are both not good. Many fans enjoy this versatile temperature & humidity meter on bed room, reptiles terrarium, garden soil, greenhouse, warehouse, basement, closet, wine cellar, humidior, and guitar room. Easy to use and install. Durable kickstand and dual magnets allow you to position it where they should be. The off/°C/°F switch gives you more freedom to choose On or Off, Celsius degree or Fahrenheit degree.



**ROOM/FREEZER
THERMOMETER DIGITAL**

Code : KMMEI- 446

Description

- LCD display shows the min/max temperature, current temperature, and battery level
- Ensure proper storage temperatures for your Food, Medicines, Chemicals, Vaccine, etc. Useful for Room, Outdoor, Freeze, Freezer, etc.,




**ROOM THERMOMETER
WOODEN TYPE**

Code : KMMEI- 447

Description

- TRADITIONAL AND ACCURATE ROOM THERMOMETER - A traditional wooden wall thermometer to measure room temperature in the home which is accurate to +/- 1C
- CLEAN & SIMPLE DESIGN! - The classic wood wall thermometer has a dark blue temperature scale on a wooden surround.
- FANTASTIC VALUE! - When it comes to accurate room thermometers, we only offer the best quality and accuracy and this item is just that!





MEDICAL EQUIPMENT & ACCESSORIES

OXYGEN MASK FOR ADULT AND CHILD, 2.0M & 5.0M

Code : KMMEI- 451

Description

- Adjustable elastic rope with mask, convenient to get comfortable fit
- Medical grade PVC material
- Each unit comes with germ-free package; could be disinfected with ethylene oxide



NASAL CANNULA FOR ADULT AND CHILD, 2.0M & 5.0M

Code : KMMEI- 452

Description

- Length: 2.0 meters
- Medical grade PVC material
- Each unit comes with germ-free package
- Could be disinfected with ethylene oxide



NEBULIZER KIT FOR ADULT AND CHILD

Code : KMMEI- 453

Description

- Soft, clear aerosol mask with anatomical form is ideal for long term use
- Gently rolled, feathered edges with nose bridges for extra comfort.
- Swivel connector for convenience of nebulization in horizontal and vertical position ensuring patient comfort.
- Larger surface area provided by unique convex cone design, ensures maximum capillary action and eliminates medicinal wastage.
- Greater nebulization rate, nebulizes 3cc medicine within 10 minutes. Specially designed for Nebulizing the patient orally.



FOLEY CATHETER

Code : KMMEI- 460

Description

- 2-Way 100% Silicone Sterilized Foley Catheter,
- Type: Two-Way Foley
- Tip: Round with Two Opposing Eyes
- Lubrication: No, Material: Silicone (Latex-Free)
- Size: 8, 10, 12, 14, 16, 18, & 20 Ch/Fr, Length: 16",
- Balloon Capacity: 5~30ml.



POSTPARTUM HEMOSTATIC BALLON

Code : KMMEI- 468

Description

- Less contraindication.
- Less-side effect and traumatic.
- Good definite effect in decreasing postpartum hemorrhage.
- Avoid of hysterectomy
- 100% Silicone
- Sterilized individual kit packed with 50ml syringe.



Items	Size	Ballon Capacity	Length (mm)
001-06	6Ch/Fr	4-5ml	300
001-08	8Ch/Fr	3-5ml	300
001-10	10Ch/Fr	3-5ml	300
001-12	12Ch/Fr	5-10ml	400
001-14	14Ch/Fr	5-15ml	400
001-16	16Ch/Fr	5-15ml/30ml	400
001-18	18Ch/Fr	5-15ml/30ml	400
001-20	20Ch/Fr	5-15ml/30ml	400
001-22	22Ch/Fr	5-15ml/30ml	400
001-24	24Ch/Fr	5-15ml/30ml	400
001-26	26Ch/Fr	5-15ml/30ml	400

Items	Size	Ballon Capacity	Length (mm)
011-24	V24Fr	500-800ml	550

