



MEDICINE SIDE CUPBOARD

BED SIDE LOCKER (ABS)

Code : KMMEI- 151

Description

- 480x480x760/800mm (LXWXH)
- Made of ABS plastics, Light, steadiness, no edging.
- Built in towel shelf and one oddment rack.



BED SIDE LOCKER

Code : KMMEI- 152

Description

- 16"x14"x60x30"mnn (LXWXH)
- M.S Tubular frame with MS cup-board.
- Stainless steel Top.
- Epoxy Powder Coated.
- Note: Also available in Complete Stainless Steel,



BED SIDE LOCKER

Supreme Code : KMMEI- 153

Description

- Overall Size 16"L x 16" W x 32" H
- CICA sheet constructions. Having one Box and one drawer.
- 5.5. Top with three sides raise edges
- Epoxy Powder Coated.



BED SIDE LOCKER

Deluxe Code : KMMEI- 154

Description

- Overall size: 16"x14"x30" (LXWXH)
- M.S Tubular frame, M.S. drawer with MS cup-board.
- Stainless steel Top.
- AM Epoxy Powder Coated.



Side Towel Railing (SS)