

**OT TABLE FULL ELECTRIC****Code : KMMEI- 131****Description**

- Power supply 220V, 50 - 60Hz.
- Corded Hand Remote.
- Five Section Radio Translucent Top Suitable for C-ARM and X-ray Procedures.
- Five section Removable Mattress.
- Matt finished polished in high quality Stainless steel Base cover, Column cover and Top frame.
- Lateral Guide rails for attachment of accessories in Stainless steel.
- Detachable and Interchangeable Head and Leg Section.
- Floor locking provided on base of the table.
- Powered Height Adjustment from 770mm - 1020mm
- Powered Trendlenburg / Reverse Trendlenburg from 20° - 20°
- Powered Lateral Tilt Right / Left from 24° - 24°
- Powered Flex and reflex position from 80° - 15°
- Manually Adjustment of Leg section to 90° down. Split away leg section to 180° as optional accessory
- Manually Adjustment of Head Section to 90° down - 90° up
- Manually adjustment by crank of Kidney Bridge to 135mm
- Standard Accessories: Side support, shoulder support, L frame, lithotomy crutches and handrest.
- Overall size: 1990 L x 575 W mm



## OT TABLES &amp; DELIVERY TABLES &amp; EXAMINATION TABLES

## OT TABLE

Code : KMMEI- 132 - SH



## Description

- Five Section SS Top with Removable Mattress.
- Stainless steel Base cover and Column cover.
- Lateral Guide rails for attachment of accessories in Stainless steel.
- Detachable and Interchangeable Head and Leg Section.
- Floor locking provided on base of the table.
- Height Adjustment by foot operated hydraulic Pump from 740 - 965mm.
- Manually adjustment of Trendlenburg / Reverse Trendlenburg from 22° - 24°
- Manually adjustment of Lateral Tilt Right / Left from 24°
- Manually adjustment of Flex and reflex position from 80° - 15°
- Manually adjustment of Leg section to 90° down.
- Manually adjustment of Head Section to 90° up and 90° down
- Manually adjustment by crank of Kidney Bridge to 100mm.
- Standard Accessories: Side support, shoulder support, L frame, lithotomy crutches and handrest.
- Overall size: 1950 L x 500 W mm

Code : KMMEI- 132-EL



**OT TABLES & DELIVERY TABLES & EXAMINATION TABLES**

**OT TABLE**

**Code : KMMEI- 133**

**Description**

- Five Section Radio Translucent Top Suitable for C-ARM and X-ray Procedures.
- Five section Removable Mattress.
- High quality Stainless steel Base cover, Column cover and Top frame.
- Lateral Guide rails for attachment of accessories in Stainless steel.
- Detachable and Interchangeable Head and Leg Section.
- floor locking provided on base of the table
- Height Adjustment by foot operated hydraulic Pump from 740 - 965mm.
- Manually adjustment of Trendlenburg / Reverse Trendlenburg from 22° - 24°
- Manually adjustment of Lateral Tilt Right / Left from 24°
- Manually adjustment of Flex and reflex position from 80° - 15°
- Manually adjustment of Leg section to 90° down.
- Manually adjustment of Head Section to 90° up and 90° down
- Manually adjustment by crank of Kidney Bridge to 135mm.
- Standard Accessories: Side support, shoulder support, L frame, lithotomy crutches and handrest.
- Overall size: 1990 L x 575 W mm



Floor Type



**OT TABLE ORTHO ATTACHMENT**

**Code : KMMEI- 133-06**

- Orthopedic Section
- Pelvic Rest
- Pelvic Support (Radio Translucent)
- Inner thigh support (Radio Translucent)
- Foot Plate with Ortho shoes
- Tibia Attachment
- Hand Surgery table (Radio Translucent)



Hanging Type



Hand Surgery



Spinal Attachment (Optional)

**OT TABLES & DELIVERY TABLES & EXAMINATION TABLES**

**OPERATION TABLE (Hi-Lo)**

**Code : KMMEI- 134**

**Description**

- M.S. Tubular Frame work on 10 cms castors.
- Height adjustment by crank mechanism with lever
- Five section SS tops
- All positions viz, 'straight', trendlenburg & reverse trendlenburg positions, Kidney bridge, Lithotomy position.
- Overall size 72" X 20" X 32"
- Fine Finished



**OPERATION TABLE IMPROVED MODEL**

(Fixed Height)

**Code : KMMEI- 135**

**Description**

- M.S. Tubular Frame work on 10 cms castors.
- Height adjustment by crank mechanism with lever
- Five section SS tops
- All positions viz, 'straight', trendlenburg & reverse trendlenburg positions, Kidney bridge, Lithotomy position.
- Overall size 72" X 20" X 32"
- Fine Finished

**LABOUR TABLE IMPROVED MODEL**

(Fixed Height)

**Code : KMMEI- 136**

**Description**

- Approx. Dimensions: 72"L x 27"W x 30"H.
- Tubular frame work mounted on PVC stumps.
- Stainless Steel top in three sections.  
U-cut at foot end.
- Trendlenburg positions adjustable on screw mechanism by handle.
- Backrest section adjustable on ratchet system.
- A pair of Lithotomy Rods with straps with height adjustable.
- Complete with self locking traction handle at side & I.V. Rod.
- Fine Finished



**LABOUR / DELIVERY TABLE TELESCOPIC**

(Fixed Height)

**Code : KMMEI- 137**

**Description**

- Approx. Dimensions: 72"L x 27"W x 30"H.
- Tubular frame work mounted on PVC stumps.
- Stainless Steel top in three sections.
- Trendlenburg positions adjustable by mechanical gear system with lever handle.
- Backrest section adjustable on ratchet .
- Leg end section can slide under the main section.
- A pair of SS Lithotomy Rods with straps with height adjustable.
- With Stainless Steel Sink/Bowl.
- Fine Finished.





**OT TABLES & DELIVERY TABLES & EXAMINATION TABLES**

**S.S. DELIVERY TABLE**

**Code : KMMEI- 138**

**Description**

- Complete stainless steel 40mm square tube frame minimum 2.0mm thickens
- Weight capacity : 150 kg - 200 kg
- Mattress 48 density PU foam 50mm thick with vinyl cover
- Three section table, Knock down, cartoon packing with complete accessories.



**OBSTETRIC DELIVERY TABLE TELESCOPIC**

(Adjustable Height) Hydraulic

**Code : KMMEI- 139**

**Description**

- Approx. Dimensions: 72"1. x 27"W x 24"-32"H.
- Tubular frame work mounted on swiveling castors.
- Stainless Steel top in three sections.
- Trendlenburg positions adjustable by pneumatic gas spring system.
- Hydraulically Height adjustable by foot paddle from side.
- Backrest section adjustable on ratchet .
- Leg end section can slide under the main section.
- A pair of knee crutches with height adjustable.
- Epoxy Powder coated.



**DELIVERY TABLE**

**Code : KMMEI- 140**

**Description**

- Approx. Dimensions: 1680L x 600W x 800H mm
- Material MS 1.2mm thickness pipe, frame made with epoxy coated.
- Pink color mattress
- Knock down.
- Cartoon Packing per piece.



**DELIVERY TABLE FULL S.S.**

**Code : KMMEI- 140 SS**

**Description**

- Stainless Stee Women Examination Bed
- Size : 1860 x 500 x 770 mm

**OT TABLES & DELIVERY TABLES & EXAMINATION TABLES**

**EXAMINATION TABLE**

**Code : KMMEI- 141**

**Description**

- Back section adjustable to 75 degree
- Powder coated steel frame, firm and anti-rust
- With a shelf to assemble non woven, easy to clean and more anti-infection
- High quantity leather, it is waterproof
- Size: 190 X 80 X 45cm
- Other size and material offered



**GYANE EXAMINATION TABLE (DELUXE)**

**Code : KMMEI- 142**

**Description**

- Frame is made up of CRCA rectangular pipe.
- Back rest adjustable by ratchet mechanism.
- Provided with parineal cut
- All metal components are pretreated and epoxy powder coated.
- Dimension : 1830 L X 560 W X 810 H mm (Approx.)
- Optional : Available with two section suitable mattress.

**EXAMINATION COUCH**

**Code : KMMEI- 143**

**Description**

- Made of CRCA sheets.
- Back rest adjustable with pneumatic pump.
- Padded with foam top in high quality washable Rexene.
- Provided with three drawers and three utility cabinet.
- Provided with foldable BP apparatus tray on the side.
- Provided with inbuilt foot step.
- Provision of saline stand.
- The complete metal component is pre treated and epoxy powder coated.



**EXAMINATION TABLE PLAIN**

**Code : KMMEI- 144**

**Description**

- Size : 72" x 20" x 30"(LxWxH).
- 2 Sections. Head rest adjustable.
- M.S. Tubular legs with rubber Stumps.
- Epoxy Powder Coated.

**GYNAE EXAMINATION TABLE**

**Code : KMMEI- 145**

**Description**

- Size : 72" x 24" x 30"(LxWxH).
- M.S. Top in 2 sections. Head rest adjustable.
- Perineal-U-cut a foot end with lithotomy rods & straps
- M.S. Tubular legs with rubber Stumps.
- Epoxy Powder Coated.

