

**OXYGEN CONCENTRATOR**



**OXYGEN CONCENTRATOR 5L**

Code : KMMEI- 255

**Description**

- Flow Specifications Total 0LPM - 5LPM
- Oxygen Concentration 93% ± 3% (5L/min)
- Input Power 350 VA
- Average Noise Level ≤50dB(A)
- Dimensions (mm) 340\*300\*650mm
- Weight 19 KG
- Warranty 2years



- Unique nebulization nozzle design
- Flexible adjustment of flowmeter, accurate flow data
- Easy to replace humidifier bottle.

3 kinds of alarm light.  
(Low concentration, low pressure, power outage)

LCD intelligent display



**OXYGEN CONCENTRATOR 10L**

Code : KMMEI- 256

**Description**

**Oxygen Concentration**

- Purity of oxygen reaches 93.96±3.96

**Multiple Ways to use**

- It comes with two flow-meters having flow from 0-5 liter and 0-10 liter per minute, it can deliver 5LPM oxygen with two persons at the same time, and can deliver 10LPM with one person
- It also has atomization function

**Safety**

- The machine has audiovival alarm for power failure, low oxygen level, low pressure
- Dimensions 390\*350\*720mm
- Weight 29.5KG
- Flow Specifications 10LPM (Dual 5 LPM each)
- Oxygen Concentration 93%±3%
- Power Consumption 750 Watts
- Sound Level 60dB(A)



10L Flow-meter 5L Flow-meter



Large LCD display