

STRETCHER AND WAITING CHAIRS

3 SEATER WAITING CHAIR, PU TOP

Code : KMMEI- 198

Description

- Size : 1800 x 630 x 800 mm
- Material : Stainless Steel finish
- PU on top



3 SEATER WAITING CHAIR, SS FINISH

Code : KMMEI- 199

Description

- Seating Capacity : 3 Seater
- Material : Stainless Steel finish
- Size 1500 x 500 x 900 mm
- Material: MS with SS finish



CPR BOARD

Code : KMMEI- 200

Description

- The CPR Board is ideal for CPR training and CPR administration by providing a firm surface that helps to properly position a patient.
- It also helps to ensure a constant open airway. The head cup maintains proper alignment of a patient's head the contoured sides make lifting the patient easier.
- Dimensions 60 x 45 x 18 cm; 1.3 Kilograms (HDPE)



SPINAL BOARD

Code : KMMEI- 201

Description

- Size : 184 x 45 x 7 cm (L x W x H)
- Self Weight : 7.5 kg
- Load Weight : 159 kg



FOLDING STRETCHER SINGLE FOLD ALUMINIUM

Code : KMMEI- 202 AL

Description

- Product size : 210 x 53 cm.
- Frame made of high quality aluminium.
- The fabric is 600D polyester covered PVC.
- Weight capacity 160kg.



STRETCHER AND WAITING CHAIRS

FOLDING STRETCHER SINGLE FOLD MILD STEEL

Code : KMMEI- 202 MS

Description

- Product size : 210 x 53 cm.
- Made of oxford fabric covered PVC, waterproof.
- Weight capacity 160kg.
- Option : Also available in PVC.



FOLDING STRETCHER DOUBLE FOLD ALUMINIUM

Code : KMMEI- 203 AL

Description

- Product size : 210 x 53 cm.
- Made of oxford fabric covered PVC, waterproof.
- Weight capacity 160kg.
- Option : Also available in PVC.



FOLDING STRETCHER DOUBLE FOLD MILD STEEL

Code : KMMEI- 203 SS

Description

- Product size : 210 x 53 cm.
- Made of oxford fabric covered PVC, waterproof.
- Weight capacity 160kg.
- Option : Also available in PVC.



DEAD BODY BAG

Code : KMMEI- 204

Description

- Product size : 220x90 cm.
- Made of oxford fabric covered PVC, waterproof.
- Weight capacity 160kg.
- Option : Also available in PVC.