

SURGICAL INSTRUMENT

BASIC PHYSIOLOGY KIT

Code : KMMEI- 475



EPISIOTOMI KIT

Code : KMMEI-476



DENTIST DAILY KIT

Code : KMMEI-477



DIAGNOSTIC SET

Code : KMMEI- 478



DRESSING KIT

Code : KMMEI- 479



HOME BANDAGE KIT

Code : KMMEI- 480




FIRST AID POUCH

Code : KMMEI- 481



SURGICAL KIT


Code : KMMEI- 482



DELIVERY

INSTRUMENT SET

Code : KMMEI- 483





SURGICAL INSTRUMENT


SUTURE KIT

Code : KMMEI- 484



VAGINAL EXAMINATINO KIT

Code : KMMEI- 485



PODIATRY STUDENTS KIT

Code : KMMEI- 486



POST-MORTEM KIT

Code : KMMEI- 487



DISPOSABLE VAGINAL SPECULUM

Code : KMMEI- 488

- Sterilized Individual Packing,
- Available in Small, Medium and Large Size



CESAREAN SECTION

Surgical Instrument Set

Code : KMMEI- 489



ADVANCE PSHYCOLOGY KIT

Code : KMMEI- 490



CIRCUMCISION KIT

Code : KMMEI- 491



ANATOMY KIT

Code : KMMEI- 492



LARYNGOSCOPE

Code : KMMEI- 493

Description

- A laryngoscope is an instrument used to view the larynx (voice box), the opening to the trachea and the lungs. It consists of a blade, which goes into the patient's mouth, and a handle. This laryngoscope set contains four differently sized Macintosh blades made of stainless steel.
- The blades of the laryngoscope are made of 304-grade stainless steel. The material makes the blade safe to use and durable
- The handle is made of durable stainless steel. It sports tiny grooves all over its body to allow a good grip
- An LED lamp is fitted in the blade which helps you view the throat clearly. The lamp is detachable and replaceable
- There are different sizes of blades that can be detached and attached to the handle according to your professional requirement.



SURGICAL INSTRUMENT

VIDEO LARYNGOSCOPE

Code : KMMEI- 496

Description

CMOS CAMERA WITH LED - 2MP CMOS camera, High Definition Anti Fog lens and integrated LED light source for high performance visualisation DISPOSABLE BLADES - Affordable, sterile, vacuum packed disposable blades connects securely to the video laryngoscope handle with a lock. 4" INCH LCD DISPLAY- High resolution 4' anti-Reflective LCD display designed for maximum visibility ERGONOMIC HANDLE- Anti microbial, ergonomically designed reusable handle to assist clinicians while use 180° SCREEN ROTATION- 180° Horizontal Screen Rotation along with Vertical movements for enhanced clinical assistance QUICK DATA TRANSFER - Videos and images can be transferred via a micro USB port for easy storage and access BATTERY INDICATOR - Rechargeable Lithium ion battery with a 2Hr capacity with a low battery indicator BUILT IN MEMORY - 4GB built in memory, for continuous video recording and image capturing with an use of a button.



PLASTER CUTTER ELECTRIC DRILL MACHINE TYPE

Code : KMMEI- 497

Description

- Supporting To the Doctor's Hand and Patient's Limbs & Soft Tissues. Suitable for cutting POP & Synthetic casting tape
- Oscillating type Cutting Action; So it will never harmful if blade touch the patient's tissue.
- Designed & manufactured for use in OPD & Casualty. Specially for Medical & Surgical Professionals. - The blade oscillates very slightly thousands of times per minute
- Easy to remove the plaster bandages POP & Synthetic casting tape. Smoothly Works on Regular AC Current & Fuse given for extra safety. Easily fix and remove blade.
- Drill Machine Type,



PLASTER CUTTER ELECTRIC METAL BODY

Code : KMMEI- 498

Description

- Supporting To the Doctor's Hand and Patient's Limbs & Soft Tissues. Suitable for cutting POP & Synthetic casting tape
- Oscillating type Cutting Action; So it will never harmful if blade touch the patient's tissue.
- Designed & manufactured for use in OPD & Casualty. Specially for Medical & Surgical Professionals. - The blade oscillates very slightly thousands of times per minute
- Easy to remove the plaster bandages POP & Synthetic casting tape. Smoothly Works on Regular AC Current & Fuse given for extra safety. Easily fix and remove blade.
- Long Body Grinding Machine Type,

