

**TROLLEY**

**INSTRUMENT TROLLEY ALL S.S**

**Code : KMMEI- 155**

**Description**

- Made up of SS pipes and sheets mounted on 10 mm castors.
- Two SS shelves, top self with 3 side railing and lower self with 4 side railing.
- Dimension: A) 600L x 460W x 810H mm  
B) 760L x 460W x 810H mm  
C) 910L x 460W x 810H mm



**INSTRUMENT TROLLEY**

**Code : KMMEI- 156**

**Description**

- Overall size : 60L X 45W X 80cms.H.
- M.S. Tubular construction with SS top & shelf.
- Mounted on 50 mm dia. Castors.
- Three side guard railing on upper shelf.
- Epoxy Powder Coated



**INSTRUMENT TROLLEY**

(Knock down construction)

**Code : KMMEI- 157**

**Description**

- Overall size : 76L X 45cris.W X 81cms.H.
- Stainless Steel Tubular construction.
- Stainless Steel Top & shelf.
- Four side guard railing on both shelves.
- Mounted on 100mm dia. Castors
- With Stainless Steel bucket & basin & push handle.



**DRESSING TROLLEY ALL S.S.**

**Code : KMMEI- 158**

**Description**

- Overall size : 76L X 45cris.W X 81cms.H.
- Stainless Steel Tubular construction.
- Stainless Steel Top & shelf.
- Four side guard railing on both shelves.
- Mounted on 100mm dia. Castors
- With Stainless Steel bucket & basin & push handle.





**TROLLEY**

**DRESSING TROLLEY**

**Code : KMMEI- 159**

**Description**

- Overall size 76L X 45cms X 81cms
- M.S. Tubular construction.
- Stainless Steel Top & shelf
- Four side guard railing on both shelves.
- Mounted on 100mm dia. Castors,
- With Stainless Steel bucket, basin & push handle.
- Epoxy Powder Coated



**MAYO'S TROLLEY ALL S.S.**

**Code : KMMEI- 160**

**Description**

- Overall size : 76L X 46W X 81H cms.
- Tray size 5.5. : Appox.:56Lx 40W cms
- Frame work S.S. Rectangular tubes mounted on 5 cms castors, height adjustable by side Knob.

Note: Also available in M.S. frame with SS tray.



**MAYO'S TROLLEY WITH GEAR**

**Code : KMMEI- 161**

**Description**

- Made up CRCA pipes mounted on 5 cms castors.
- Height adjustable by gear provided on side.
- Dimension of SS tray : 600 L x 460W mm.



**MEDICINE TROLLEY WITH FOUR DRAWERS**

**Code : KMMEI- 162**

**Description**

- Size 28" x 18"W x 32" H.
- M.S. Tubular frame mounted on 10 cms dia. Castors
- Stainless Steel top & Shelf.
- Two Drawers under each shelf.
- Three side guard rail on upper shelf.
- Epoxy Powder Coated



**TROLLEY**

**MEDICINE TROLLEY WITH TWO DRAWERS**

**Code : KMMEI- 163**

**Description**

- Size : 650Lx460Wx900H mm
- S.S. Tubular frame mounted on 10mm dia. castors. Stainless Steel top & shelf.
- Two drawers under top shelf.
- Four side guardrail on both shelves
- Cartoon Packing



**MEDICINE SIDE CUPBOARD**

**Code : KMMEI- 164**

**Description**

- Made up of CRCA pipe/sheets with SS top.
- Provided with two shelves and two drawers.
- Mounted on 5 cm casters.
- Pretreated and epoxy powder coated.
- Approx. Dimension: 405mm
- L\* 405 mm W\* 810 mm H.



**STRETCHER TROLLEY ALL S.S.**

**Code : KMMEI- 165**

**Description**

- Overall Size: 213 x 56 x 381 cms (LXWXH)
- Stainless Steel Tubular construction with SS removable stretcher.
- Mounted on 150mm dia. Coaster Two with brakes)
- Both side IV attachment with one SS I.V, Double Hook Rod.



**STRETCHER TROLLEY (DELUXE)**

**Code : KMMEI- 166**

**Description**

- Removable with Fixed Height
- Frame is made of CRCA tube.
- Trolley mounted over 150mm caster for ease of movement.
- Provision of saline stand.
- Provided with swing down railing.
- Provided with cylinder cage
- \*Dimension (Approx) : 2100 L x 560 W x 810 H



**STRETCHER TROLLEY (STANDARD)**

**Code : KMMEI- 167**

**Description**

- Removable with Fixed Height
- Frame is made of CRCA tube.
- Trolley mounted over 150mm caster for ease of movement.
- Provision of saline stand.
- Dimension (Approx) 2100 L x 560 W x 810 H
- Note : Optional mattress





**TROLLEY**

**CRASH CART TROLLEY STAINLESS STEEL**

**Code : KMMEI- 168**

**Description**

- Overall approx. Size: 940mm Lx 490mm W x 1535mm H
- Stainless Steel Tubular frame.
- Six removable plastic bins and two polystyrene lockable storage Units with three drawers each.
- Four swivel castors, 150mm dia. (Two with brakes)
- Complete with corner buffers, oxygen cylinder holder, S.S IV Rod, S.S Top & shelves, cardiac massage board and electric lamp.



**CRASH CART TROLLEY**

**Code : KMMEI- 169**

**Description**

- Overall approx. Size: 940 mm L x 490 mm W x 1535 mm H
- M.S. Tubular frame.
- Six removable plastic bins and two polystyrene lockable storage units with three drawers each.
- Four swivel castors, 150mm dia. (Two with brakes)
- Complete with corner buffers, powder coated oxygen cylinder holder, S.S IV Rod, S.S. Top & shelves, cardiac massage board and electric lamp.
- Epoxy Powder Coated



**ATTENDANT CHAIR BED**

**Code : KMMEI- 170**

**Description**

- Size : 1900Lx620Wx580H mm
- Made of powder Coated mild Steel.
- P.U. foam top.
- Weight capacity 160kg
- Foldable.

