



WHEEL CHAIR, WALKER & WALKING STICKS

WHEEL CHAIR

Code : KMMEI- 241

Description

- M.S. frame chrome finish.
- Detachable foot rest.
- Aluminium foot plate.
- Foldable back rest.
- Weight capacity 100 - 125kg.



WHEEL CHAIR WITH UNITED BREAK

Code : KMMEI- 241 B

Description

- M.S. frame chrome finish.
- Detachable foot rest.
- Aluminium foot plate.
- Foldable back rest.
- Weight capacity 100 - 125kg.



WHEEL CHAIR FOR CHILD

Code : KMMEI- 242

Description

- Standard children wheel chair.
- M.S. frame green powder coating.
- Fixed foot rest
- Black plastic foot plate.
- Weight capacity 75kg



WHEEL CHAIR WITH COMMODE

Code : KMMEI- 243 E

Description

- M.S. frame chrome finish.
- Aluminium foot plate.
- Weight capacity 100 - 125kg.



WHEEL CHAIR, WALKER & WALKING STICKS

WHEEL CHAIR WITH CAMMODE

Code : KMMEI- 243-R

Description

- Reclining back rest 900-1800.
- M.S. frame chrome finish.
- Detachable arm rest.
- Detachable space and elevating foot rest,
- Aluminium foot plate.
- Weight capacity 100 - 125 kg.



WHEEL CHAIR (Aluminium)

Code : KMMEI- 244

Description

- Aluminium wheel chair with foldable back rest.
- Aluminium frame with silver liquid coating.
- Black plastic foot plate.
- Weight capacity 100 - 125kg.
- Brake : United Brakes.



THREE IN ONE WALKER

Code : KMMEI- 245

Description

- Foldable Commode Chair, Walker, Shower Chair,
- Aluminium Frame, Color: Oxidant bright Silver,
- Size: 5101 x 430W x 470-650H (mm),
- Weight Capacity 75 -90 KG



COMMODOE CHAIR

Code : KMMEI- 246

Description

- Size: 520Lx550Wx450-550H (mm)
- Weight Capacity 75 KG,
- Foldable Aluminium Commode Chair
- Frame color: Oxidant bright silver.



CAMMODE CHAIR DELUXE MODEL

Code : KMMEI- 246 B

Description

- Weight Capacity: 75 KG,
- Cushioned Seat and Back Support for Comfort,
- Foldable Chair, with Removal Pan with Lid

SIZE	62	50	45	93	42	52	75	80*27*59	2	7	508
MODEL											



WALKER & WALKING STICKS

ALUMINIUM FOILDING WALKER

Code : KMMEI- 247

Description

- Overall size: 490Lx 590W x Adjustable height 760-940mm
- Weight capacity : 75kg
- 4 pieces packing carton.



ADJUSTABLE WALKING STICK

Code : KMMEI- 248

Description

- Aluminium adjustable stick
- Available in different color and shape.



ALUMINIUM CRUTCHES

Code : KMMEI- 249

Description

- Adjustable Height Crutch
- Available in Small, Medium and Large Size,
- Soft Cushion, and Soft Hand Grip, Nylon Stump,
- Light Weight, Maintenance free,